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- CG.1
- CG.49
- CG.61.1
- CG.61.2
- CG.103

Unconventional prosthesis for reconstruction of upper limbs

This surgical technique was developed by

Technology Division - Baumer S.A. - Rev. 4.001

This technique covers the following products:

- Proximal humerus;
- Proximal modular humerus;
- Humeral diaphyseal;
- Total humerus and elbow articulated;
- Distal humerus and elbow articulated.

Figure 1:
306.REF
(Proximal Humerus)

Figure 2:
343.REF (Proximal
Modular Humerus)

Figure 3:
1265.REF
(Humerus Diaphyseal)

Figure 4:
1264.D/E (Total
humerus and elbow
articulated)

Figure 5:
350.D/E (Distal
humerus and elbow
articulated)

Figure 6:
1350.D/E (Joint Total
Elbow)

Figure 7:
338.REF (Original
Radio Proximal
M.H.M.)

Indications

The indication, purpose or use for which the medical product is intended is for cases of partial or total replacement or reconstruction of bone segments of the upper limbs.

Due to diseases, such as:

- Replacement of bone segment due to tumor resection surgery.
- Partial or total replacement of the bone segment due to other causes.
- Joint replacement.
- Review of primary arthroplasties.
- Bone segment reconstruction.
- Trauma with substantial bone loss

Contraindications

The contraindications below should be strongly considered by the surgeon. Such contraindications include, but do not limit, the use of the Unconventional Prosthesis for Reconstruction of Upper Limbs, in cases of:

Recent history of systemic or localized infection.

- Signs of local inflammation;
- Bone immaturity;
- Severe osteoporosis;
- Morbid obesity;
- Fever or leukocytosis;
- Sensitivity or allergy to the metal of the medical product;
- All those known for the use of this medical product;
- Patients unable to follow postoperative care and instructions.

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- The use of the Unconventional Prosthesis for Reconstruction of Upper Limbs has offered means of replacement or partial or total reconstruction of bone segments of the upper limbs. However, this medical product is designed solely to replace or reconstruct the damaged bone segment and not the normal structures of the human skeleton;
- Factors such as the patient's weight, activity

level and attention to instructions regarding post-surgical limitations have effect on load and number of cycles at what medical product is exposed;

- The surgeon must be fully aware not only of the medical and surgical aspects of the medical product, but also be aware of the mechanical and metallurgical aspects of it, for the correct use of the product and guidance of the patient.
- The correct selection of the medical product is the responsibility of the doctor.
- The Unconventional Prosthesis for Reconstruction of Upper Limbs cannot withstand levels of activity and /or loads equal to those supported in normal and healthy bones.
- Although with all the scientific and technological development, metals, and their alloys, when implanted in the human body, are subject to constant environmental changes (concentration of salts, acids, and alkaloids), which can cause damage to the implant. The use of dissimilar metals in contact with each other can accelerate the corrosion process, which in turn accentuates the fatigue fracture of the medical product.
- The mixture of different biomedical materials in the same patient, except when they are compatible or similar, can increase the risk of corrosion.
- The surgeon should be attentive when performing magnetic resonance imaging, if the patient already has implants made of stainless steel. The presence of these implants can distort the results of the examination, impairing the correct diagnosis of diseases.
- The Unconventional Prosthesis for Reconstruction of Upper Limbs can loosen, break, suffer corrosion, cause pain, or even weaken the bone, particularly in young and active patients.
- The patient should be alerted to the limitations of the medical product, and that the level of physical activity, weight and overload has implied premature failure, deformation, or breakage.
- Excessive physical activity, and trauma affecting the replaced joint have resulted in premature failure to replace the bone segment, either by loss, fracture or wear of the medical product and its components ancillary. The patient should be alerted to maintain their activities according

to their situation, protecting the replaced joint against excessive "stress".

Introduction

Defects created by tumor resection or other causes should be treated on the basis of orthopedic principles. In oncology, there is no logical justification for a poor reconstruction simply because the patient is a carrier of a tumor, even if malignant.

In general, patients with primary tumors are young and due to this, more active, which submits to an extremely critical test, the design of the prosthesis, the correct selection of the implant and the method of fixation. When we decide to use a stent, we must choose an implant that is able to promote an extremely long period of use, without mechanical complications.

The main factors that influence the method of reconstruction are:

- Relief of pain caused by the presence of tumor or pathological fracture.
- Stability that allows the patient to move to activities of daily living.
- Mobility allowing the best possible function of the joint.
- Durability of reconstruction compatible with the increasingly long survival of patients with skeletal neoplasms.
- Emotional and psychological impact that would be caused by a mutilating surgery.

Aggressive benign lesions or low-grade malignancy of the bone can often be successfully treated through local resection and reconstruction of the affected limb. Even primary lesions of a high degree of malignancy, or metastases of bone or soft tissues, can be treated with limb preservation, especially after the institution of adequate chemical or radiotherapy treatment, always without prejudice to the patient's survival.

We have developed a new concept of unconventional prostheses for each of the regions, upper limbs, and lower limbs. These stents are based on standard measurements and are previously constructed, being available in hospitals and allowing surgeries to be scheduled immediately after diagnosis and decision of the best conduct for each case.

Currently, the use of stents is not only limited to neoplastic conditions, but also offers.

a new option in the treatment of many non-neoplastic conditions, involving bone and joint loss caused by trauma, osteo-metabolic diseases, and failures of conventional arthroplasties.

Surgical Planning for the Unconventional Prosthesis for Reconstruction of Upper Limbs

To prescribe and request a stent graft does not conventional radiographs were used, in anteroposterior and lateral view, of the entire bone to be operated and often of the contralateral bone, when the bone integrity is altered by the tumor to the point that a reliable measurement is difficult.

With the use of "templates" (transparency – when applicable), by the method of overlapping in radiographic examination, it is possible to predetermine the size and model of the Unconventional Prosthesis for Reconstruction of Upper Limbs to be utilized.

Attention: Ask the Radiology Service of your Hospital for the full-size X-ray(s) (1:1 scale). Radiographs in scales other than 1:1 or reduced digital may compromise planning.

Note: To obtain the Magnification Factor of 1:1 - The distance between the X-ray source and the image (DFI) must be equal to 100 cm. The radiopaque millimeter (or centimeter) ruler should be next to the bone, in the longitudinal direction, and should not be supported on the examination table, which would enlarge its image when compared to the bone.

In this radiography, the measurements studied are:

- The length of the area to be respected, which corresponds to the length of the body of the endoprosthesis.
- The length of the remaining normal diaphysis, to evaluate the maximum length of the nail.
- The diameter of the remaining normal medullary canal, which determines the diameter of the stent rod.
- The type of joint used

Access Route

The non-conventional prosthesis can be implanted by the access route of the surgeon's choice.

Resection

The length of the area to be resected will correspond to the length of the "Endoprosthesis Body" determined in the surgical planning.

The surgeon must, first, based on imaging tests, determine which region affected by the pathology he will address.

According to the pathology, the surgeon may opt for total or partial replacement of the humerus and/or elbow.

Preparation of the humerus

Once the endoprosthesis that best suits the case is defined, start the resection of the compromised segment of the humerus and then regularize the surface according to the sequence of images below:

F. Level of resection planned by the Surgeon

G. Osteotome marking of a line showing the center of the anterior surface of the humerus. This line will serve as a reference guide to guide the rotation of the prosthesis and avoid inadequate internal or external rotation.

preserving the humeral condyles, proceeding to the modeling osteotomies using osteotome or sagittal orthopedic saw, as illustrated in the image above.

H. After osteotomy, the bone end should be regularized with the top cutter.

Preparation of the ulna and radius (elbow)

Once the endoprosthesis that best suits the case is defined, start with the removal of the end of the olecranon with the aid of an oscillating saw, the amount of bone removed should allow access to the milling of the medullary canal, as shown below:

D. Removal of the end of the olecranon

I. Next, the milling of the medullary canal of the bone portion to the correct diameter of the stem of the endoprosthesis should be started, taking into account the technique of bone cementation. Start the procedure with smaller diameter milling cutters, progressively progressing with immediately larger milling cutters. Great care must be taken not to make a false path, which would cause the extrusion of the cement and the prosthesis rod.


J. For the implantation of the articulated elbow endoprosthesis (Code: 1350.D/E) in the distal region of the humerus (elbow joint) one can choose to implant the component

E. The removal of the radius head is necessary to start the milling of the medullary canal of the same, as shown above.

Note.: For the Original Proximal Radio Prosthesis M.H.H (338.REF) the bone removal must correspond to the length of resection of the prosthesis (9, 12 and 18mm). Instrumental tests are available for verification of bone resection. As shown in the image below.

W (mm)
9,0
12,0
18,0

F. Next, the milling of the medullary canal of the bone portion to the correct diameter of the stem of the endoprosthesis should be started, considering the technique of bone cementation. The procedure should always be started with smaller diameter milling cutters, progressively progressing with immediately larger milling cutters.

Instructions for assembly and disassembly of the definitive implant

The correct handling of the implants is extremely important. To ensure proper implantation

of the prosthesis, it is not recommended to disassemble the product, as it may contribute to early failure of the implant. If necessary, follow the instructions below.

Figure 8:
Refers to the Distal Humerus and Elbow Original Baumer Articulated Endoprosthesis

Figure 9:
Refers to the Distal Humerus and Elbow Original Baumer Articulated Endoprosthesis

The sequence of the technique of disassembly is presented below:


III) Dismantling procedure

To disassemble the prosthesis, perform the following operations:

- Loosen the latch screw from the shaft with the aid of a 2.5mm Allen (hexagonal) wrench.
- Release the locking screw with two counterclockwise turns (Warning: It is not necessary to remove the screw completely).

c. Components of the prosthesis:

- 1-Proximal Ulna Component
- 2-Proximal Radio Component
- 3-Shaft Lock Screw
- 4- Component distal of Humerus

The Sequence of the assembly technique is presented below


IV) For the assembly of the prosthesis perform the following operations:

d. Position the proximal component of Ulna and Radius in the distal component of Humerus in the impaction hole.

e. Make sure that the components of the Ulna and Radio are positioned correctly (Attention: The hole of the ulna component must be aligned with the latch screw).

f. With the hexagonal wrench tighten the locking screw and check the blocking of the internal and external rotation movements. Confirm the height of the screw-lock. It should be aligned to the surface of the assembly. This visual check is another way to ensure that the components are assembled correctly.

Attention: If the screw is spare as in the image above it is in the incorrect position and will not lock the joint properly.

Instruments for prostheses not conventionalis

To perform non-conventional prosthesis surgery for upper limbs it is recommended to use the Baumer instrument kits with the following codes:

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